

		GINGIVAL HEIGHT OF ABUTMENT								
		1 mm	2 mm	3 mm	4 mm	5 mm	6 mm	Ang 1,5		
COMPATIBLE IMPLANT CONNECTION with		PRODUCT CODES								
ASTRA TECH	Lilac									
MOZO GRAU	InHex WP									
NEOBIOTECH	IS 3.5/4.0/4.5/5.0									
BIOGENESIS	Aticon 4.5									
COWELLMEDI	Submerged	DS50451	DS50452	DS50453	DS50454	DS50455	DS50456			
DIO	UFII 3.8/4.0/4.5/5.0									
DENTIS	S-clean									
MEGAGEN	AnyOne Internal									
MEGAGEN	EZ Plus 4.0/4.5/5.0									
ASTRA TECH	Aqua									
MOZO GRAU	InHex RP									
BIOGENESIS	ATICON 3.5	DS50351	DS50352	DS50353	DS50354	DS50355	DS50356			
EUROTEKNIKA	Naturactis.Naturall.Natea									
MEGAGEN	EZ Plus 3.5									
OSSTEM IMPLANT	TS 3,5	DS54351	DS54352	DS54353	DS54354	DS54355	DS54356			
OSSTEM IMPLANT	TS 4.0/4.5/5.0	DS54451	DS54452	DS54453	DS54454	DS54455	DS54456			
MIS	Internal 3,3 Narrow platf.									
JDENTALCARE	Evolution S 3.2	DS20331	DS20332	DS20333	DS20334	DS20335	DS20336			
PALTOP	Internal hex. Narrow 3.25									
ZIMMER	Internal Hex. SP									
ALPHA-BIO	Internal Hex. SP									
BIO-HORIZON	Internal 3,5									
IMPLANT DIRECT	Legacy 3.7									
MIS	Internal 3.75/4.2 SP	DS20351	DS20352	DS20353	DS20354	DS20355	DS20356	DS2035AN1.5	DS2035AN3.0	
AB IMPLANTS	Internal Hex.									
ACE	Internal Hex. SP									
JDENTALCARE	JDEvolution									
AVINENT	HI SP									
ZIMMER	Internal Hex. WP									
BIO-HORIZON	Internal 4.0									
IMPLANT DIRECT	Legacy 4.7	DS20451	DS20452	DS20453	DS20454	DS20455	DS20456	DS2045AN1.5	DS2045AN3.0	
MIS	Internal 5.0 WP									
ACE	Internal Hex. WP									
AVINENT	HI WP									
NOBEL BIOCARE	Branemark Ext. Hex. 3.3 NP									
OSSTEM IMPLANT	US II 3.5	DS10331	DS10332	DS10333	DS10334	DS10335	DS10336			
ACE	External Hex. NP									
ACE	External Hex RP									
AVINENT	HE 4.1									
BIOGENESIS	3ICON									
BIOMET 3i	Osseotite HE 3.75/4.00									
BTK	Evo Klassic Ext									
COWELLMEDI	External									
DIO										
MEGAGEN	Anyone External 4.0									
NEOBIOTECH	EB External									
NOBEL BIOCARE	Branemark Ext. Hex. 4.1 RP	DS10401	DS10402	DS10403	DS10404	DS10405	DS10406	DS104AN1.5	DS104AN3.0	
	Speedy 4.0 Hex. Ext.									
OSSTEM	US II 4.1									
BTI	Hexágono Externo 4.1									
IMTEC	Universal Hex. Ext.									
INTRA-LOCK	3.75/4.00 Hex. Ext.									

LIFECORE	Restore 3.75/4.00 Hex. Ext.								
MIS IMPLANTS	External SP								
MOZO-GRAU	3.75/4.00 Hex. Ext.								
ACE	External Hex. WP								
NOBEL BIOCARE	Branemark Ext. Hex. 5.1 WP	DS10501	DS10502	DS10503	DS10504	DS10505	DS10506		
OSSTEM	US II 5.1								
KLOCNER SK2 4.2	SK2 4.2		DS40422 (2.5 mm)	DS40423 (3.5 mm)	DS40424 (4.5 mm)	DS40425 (5.5 mm)			
NOBEL BIOCARE	Replace 3,5								
IMPLANT DIRECT	Replant 3,5	DS12351	DS12352	DS12353	DS12354	DS12355	DS12356		
ACE	Tri-Cam 3,5								
NOBEL BIOCARE	Replace 4,3								
IMPLANT DIRECT	Replant 4,3	DS12431	DS12432	DS12433	DS12434	DS12435	DS12436		
ACE	Tri-Cam 4,3								
NOBEL BIOCARE	Active 3,0	DS16301	DS16302	DS16303	DS16304	DS16305	DS16306		
NOBEL BIOCARE	Active 3,5								
AVINENT	Ocean conical 3,5	DS16351	DS16352	DS16353	DS16354	DS16355	DS16356		
NOBEL BIOCARE	Active 4,3								
AVINENT	Ocean conical 4,3	DS16431	DS16432	DS16433	DS16434	DS16435	DS16436		
MULTI-UNIT	Multi-Unit abutment	DS1MU01 / DS1MU02							
STRAUMANN	ITI 4.8 Tissue Level								
LIFECORE BIOMEDICAL	Stage 1								
MEGAGEN	Exfeel Internal 3.5/4.2/4.8								
BIOGENESIS	ITICON								
ANTHOGRYR	Ossfit								
BTK	8Neck INT								
COWELLMEDI	Internal								
DENTIS	I-Clean	DS30401	DS30402	DS30403	DS30404	DS30405		DS304AN1,5	DS304AN3,0
DENTIUM	Simple Line								
EUROTEKNIKA	Aesthetica								
IMPLANT DIRECT	Swissplant								
MIS MEDICAL	Mistral								
OSSTEM IMPLANT	SS 4.8								
ZIMMER	Swisspluss 4.8								
PHIBO	TSH/BNT S2	DS36201	DS36202	DS36203	DS36204	DS36205			
	TSH/BNT S3/S4	DS36341	DS36342	DS36343	DS36344	DS36345			
	TSH/BNT S5		DS36502	DS36503					
PHIBO	TSA/TSDAV S3	DS38301	DS38302	DS38303	DS38304	DS38305			
	TSA/TSDAV S4	DS38401	DS38402	DS38403	DS38404	DS38405			
	TSA/TSDAV S5		DS38502						
BIOMET 3i	Certain 3,4	DS60341	DS60342	DS60343	DS60344	DS60345	DS60346		
	Certain 4,1	DS60411	DS60412	DS60413	DS60414	DS60415	DS60416		
MEGAGEN	AnyRidge	DS51401	DS51402	DS51403	DS51404	DS51405	DS51406		
TRI Implants	Narrow 3,3	DS24321	DS24322	DS24323	DS24324	DS24325			
DENTALTECH	Internal hex. 3.25	DS22321	DS22322	DS22323	DS22324	DS22325			

References available coming soon